


# StimConnect™

ELONGATION CABLES FOR TEMPORARY STIMULATION



# Accessories Stimulation

## StimConnect™

Disposable and reusable bipolar elongation cables for connecting myocardial electrodes and catheters to external pacemakers for temporary stimulation. The cables are available in blue and white and have a length of 2 m.

### Disposable Elongation Cables

Disposable elongation cables are sterile packed for safe and immediate usage.

| Distal Side (Female)<br>Connection Myocard-Electrode/ Catheter | Proximal Side (Male)<br>Connection Pacemaker     |
|--|--|
| Terminal box with screw side socket for naked 2 mm pins        | 2 mm shrouded pins for collets<br>Osypka adapter |
| Terminal box with collets for naked 2 mm pins                  | Open 2 mm pins                                   |
| Terminal box for touch proof 2 mm shrouded pins                | Touch proof 2 mm shrouded pins                   |
|  | Medtronic adapter                                |
|  | Biotronik adapter                                |

### Reusable Elongation Cables

Reusable extension cables can be sterilized 25 times by ethylene oxide or steam. A punch card allows the identification of the sterilization cycles.

Possible distal and proximal connection variants, freely combinable:

| Distal Side (Female)<br>Connection Myocard-Electrode/ Catheter | Proximal Side (Male)<br>Connection Pacemaker |
|--|--|
| Terminal box with screw side socket for naked 2 mm pins        | Open 2 mm pins                               |
|  | Touch proof 2 mm shrouded pins               |
|  | Medtronic adapter                            |
|  | Biotronik adapter                            |

## EasyStim™

EasyStim adapter for connecting touch proofed 2 mm pins (e.g. catheter) to collets, which allow to connect naked 2 mm pins only. Double (red and black), sterile or non-sterile packed. Packaging unit: 20 double packs



### Product Code

EasyStim sterile: EasyStim non-sterile:

EasyStim.S

EasyStim

## StimConnect – Disposable Elongation Cables



**Distal Side (Female)**  
**Connection Myocard-Electrode/ Catheter**

**Proximal Side (Male)**  
**Connection Pacemaker**



Terminal box with  
 collets for naked  
 2 mm pins  
**Code: TCO**



Terminal box for  
 touch proof 2mm  
 shrouded pins  
**Code: TSP**



Terminal-box with  
 screw side socket for  
 naked 2 mm pins  
**Code: TBOX**

Packaging unit: 10 pieces

**Product Code for Disposable Elongation Cables:**

StimConnect.



Color:  
 W = White  
 B = Blue




Code:  
 Distal Side



Code:  
 Proximal Side


Osyпка adapter  
**Code: OSY**



Naked 2 mm pins  
**Code: OP**



Touch proof 2 mm  
 shrouded pins  
**Code: SP**



Medtronic adapter  
**Code: MED**



Biotronik adapter  
**Code: BIO**



Example: StimConnect.W/TCO/OSY (corresponds to the elongation cable shown above)

# Accessories Stimulation

## StimConnect – Reusable Elongation Cables




**Distal Side (Female)**  
**Connection Myocard-Electrode/ Catheter**

**Proximal Side (Male)**  
**Connection Pacemaker**




Terminal-box with  
 screw side socket for  
 naked 2 mm pins  
**Code: TBOX**

Naked 2 mm pins  
**Code: OP**



Touch proof 2 mm  
 shrouded pins  
**Code: SP**



Medtronic adapter  
**Code: MED**



Biotronik adapter  
**Code: BIO**



Packaging unit: 5 pieces

**Product Code for Disposable Elongation Cables:**

StimConnect. [R] [ ] / TBOX / [ ]

Color:  
 W = White  
 B = Blau

Code  
 Proximal Side

Example: StimConnect.RB/TBOX/OP (corresponds to the elongation cable shown above)

Osypka Medical GmbH

Albert-Einstein-Strasse 3  
 12489 Berlin, Germany  
 Phone: +49 (30) 6392 8300  
 Fax: +49 (30) 6392 8301  
 E-Mail: mail@osypkamed.com  
 www.osypkamed.com

StimConnect and EasyStim are trademarks of Osypka Medical  
 Copyright © 2019 Osypka Medical



STOMENT190612