


- Components:
- sheath
 - dilator
 - guidewire
 - puncture needle




Callisto™

INTRODUCER KIT (FEMORAL)



ORDERING INFORMATION					
Femoral Sheath			(Packaged five per unit box)		
CATALOGUE NUMBER	DIAMETER	LENGTH (cm)	NEEDLE GAUGE (Ga)	GUIDEWIRE DIAMETER (inch)	GUIDEWIRE LENGTH (cm)
CFS051135	5F	11	18	0.035	45
CFS051138	5F	11	18	0.038	45
CFS061135	6F	11	18	0.035	45
CFS061138	6F	11	18	0.038	45
CFS061638	6F	16	18	0.038	80
CFS062338	6F	23	18	0.038	80
CFS071135	7F	11	18	0.035	45
CFS071138	7F	11	18	0.038	45
CFS071638	7F	16	18	0.038	80
CFS072338	7F	23	18	0.038	80
CFS081135	8F	11	18	0.035	45
CFS081138	8F	11	18	0.038	45
CFS081638	8F	16	18	0.038	80
CFS082338	8F	23	18	0.038	80
CFS091138	9F	11	18	0.038	45




PROVIDING A SUCCESSFUL ENTRY

CallistoTM


INTRODUCER KIT (FEMORAL)

Right choice at the beginning, better finish in the end

- Well balanced flexibility and kink-resistance
- Specially rounded dilator tip to allow easy vessel entry and atraumatic intraluminal advancement




- Smooth transition between wire, dilator and sheath to minimize the vessel injury



- High versatility valve to avoid leakage of blood and air aspiration



- Quarter turn screw lock to eliminates risk of dilator back-out and allows easy unlocking and withdrawl



- Standard color of dilator and sheath hub helping users to detect the French size easily



- choices of various length of sheath

