


# Rayband™

Radial Artery Compression Tourniquet **FDA** **CE** 0197

## Clear Visualization

Easy control of hemostasis process



## 2mL Controlled Air Release Locker

Provides precise and step-wise deflation in an easy and safe way

>>> Single Air Bag Structure Tourniquet <<<



>>> Dual Air Bag Structure Tourniquet <<<



## Product Information

### Features

- Clear visualization through transparent air bag
- Controlled air release locker
- Adjustable bandage length


### Benefits

- Easy control of hemostasis process
- Convenient and safe deflation
- Comfortable hemostasis for patients

### Ranges

- Air bag structure: single / dual
- Max. air volume: 15ml / 20ml
- Connector: Luer / Push-pull
- Bandage length: 21cm/24.5cm/27cm/29cm

## Direction for Use



## Ordering Information

Catalogue No.	Air Bag Structure	Bandage Length (cm)	Max Injection Air Volume (mL)	Connector Type
LPYP15S	1	21	15	Push-pull
LPYP15	1	24.5	15	Push-pull
LPYP15L	1	27	15	Push-pull
LPYP15XL	1	29	15	Push-pull
LPYP20S	2	21	20	Push-pull
LPYP20	2	24.5	20	Push-pull
LPYP20L	2	27	20	Push-pull
LPYP20XL	2	29	20	Push-pull
YPR15S	1	21	15	Luer
YPR15	1	24.5	15	Luer
YPR15L	1	27	15	Luer
YPR15XL	1	29	15	Luer
YPR20S	2	21	20	Luer
YPR20	2	24.5	20	Luer
YPR20L	2	27	20	Luer
YPR20XL	2	29	20	Luer